










Medimec® International srl
 Via Proventa, 52 - 48018 Faenza (RA)
 Tel. 0546/46870 - Fax 0546/46467
 P.IVA 02259030399
medimec@medimec.it | www.medimec.it

Deambulatori Meyland-Smith: tabella riepilogativa componenti

	Miniwalk Mis. 1	Miniwalk Mis. 2	Meywalk 2000 Small	Meywalk 2000 Medium	Meywalk 2000 Large	Meywalk 4 Medium	Meywalk 4 Large
							
Altezza sella (senza carico)	31 - 49 cm	31 - 64 cm	48 - 66 cm	62 - 86 cm	76 - 100 cm	62 - 87 cm	83 - 106 cm
Portata massima	40 kg	40 kg	100 kg	100 kg	100 kg	130 kg	130 kg
Molla (colore)							
da 10 a 15 kg (verde)	X	X	-	-	-	-	-
da 15 a 20 kg (argento)	X	X	-	-	-	-	-
da 20 a 45 kg (nero)	X	X	-	-	-	-	-
da 20 a 45 kg (bianco)	-	-	X	X	X	X	X
da 45 a 65 kg (giallo)	-	-	X	X	X	X	X
da 65 a 80 kg (blu)	-	-	X	X	X	X	X
da 80 kg in su (rosso)	-	-	X	X	X	X	X
Anello toracico (circ. interna)							
70 cm	X	X	X	-	-	X	X
90 cm	X	X	X	X	X	X	X
105 cm	-	-	X	X	X	X	X
120 cm	-	-	-	X	X	X	X
130 cm	-	-	-	X	X	X	X