



# Woven Air

# Sprint

## Technical Data Sheet



<b>Foam Base Composition</b>	<b>EVA Closed Cell Foam</b>
<b>Foam Pre-ischial Shelf Composition</b>	Enduro Open Cell High-resilient Foam
<b>Woven Air Tube Quantity</b>	8 per cushion (10 for cushions 18x22 – 20x20)
<b>Cushion Weight</b>	900gms (18x18)
<b>Fire-retardant Fabrics &amp; Materials</b>	yes
<b>Incontinence Cover</b>	Not included. This is a separate Accessory
<b>Contour Profile</b>	20mm
<b>Cushion Height</b>	60mm (from base to under thighs)
<b>Outer Cover Composition</b>	Upper panel – FLO2™ 4-way stretch Water-resistant Fabric. Side Borders – 3d Spacer Fabric Under panel – Rubberised non-slip Urethane Fabric Includes hook receptive loop panels & a carry handle.
<b>Pelvic Obliquity Pads?</b>	Not included. This is a separate Accessory